

Food and environmental analytics

114765 Copper Test Method: colorimetric (refill pack see art. 1.14767.0001) 0.3 - 0.6 - 1.0 - 1.5 - 2 - 3 - 5 - 7 - 10 mg/l 100 Tests Microquant®

Formula Hill

Chemical formula

Molar mass

Density

pH value

Form

Color

Odeur

Melting point

Boiling point

Vapor pressure

Solubility in water

Explosion limits

Solubility in ethanol

Flash point

Solubility in chloroform

Thermal decomposition

Water absorbtion

Heat of evaporation

Saturation conc. (air)

Evaporation number

Dielectric constant

Ignition temperature

Dipole moment

Bulk density

Viscosity dynamical

Refractive index

Viscosity kinematical